Flexible Installations of Controlled Atmosphere Units

« Janny MT preservation »

Broccoli


- Storing during production peaks
- Maintaining flavour quality
- Maintaining colour and firmness
- Avoiding weight loss (less than 1% after 45 d)
- Minimize losses during long storage

Picking recommendations

- Match maturity to the desired preservation time
- Pick extremely healthy products
- Do not store products after rain
- Favour small calibres

Key figures

- 1. Use Janny MT module 610 6
- 2. Fill module with **150 kg loose**
- 3. Put in cold storage at **+0°C**
- 4. Wait until the products are cooled to heart
- 5. Close the lid
- 6. Check gas stabilisation
- 7. Possible preservation up to 45 days
- 8. Open modules
- Let the module 12 hours with open lid in cold room for rehabilitation




6 Open membranes
0 Cap on

Temperature	O ₂ rate	COarate	Preservation

	Temperature	O ₂ rate	CO ₂ rate	Preservation time
Classic Cold	+0°C	21 %	0 %	21 d
Janny MT CA	+0°C	1 to 5 %	1 to 5 %	45 d

Ideal curve

Gas content monitoring for Aspargus





Levels are read using the Tiempo Test measuring instrument



Broccoli

. . .

Sté JANNY MT SAS
La Condemine 71260 Péronne France
Tel: 0033.3.85.23.96.20
Fax: 0033.3.85.36.96.58 / 0033.3.85.23.96.29

Mail: jannymt@jannymt.com Web: www.jannymt.com