



JANNY MT



Flexible Installations of Controlled Atmosphere Units

Cardoon



Cardoon


- Storing during production peaks
- Maintaining flavour quality
- Maintaining colour and firmness
- Avoiding weight loss (less than 1% after 45 d)
- Minimize losses during long storage

Picking recommendations

- Match maturity to the desired preservation time
- Pick extremely healthy products
- Do not store products after rain

Key figures

1. Use Janny MT module 610 - 6
2. Fill module with **150 kg loose**
3. Put in cold storage at **+2°C**
4. **Wait until the products are cooled to heart**
5. Close the lid
6. Check gas stabilisation
7. Possible preservation up to **45 days**
8. Open modules
9. Let the module **12 hours with open lid in cold room for rehabilitation**




6 Open membranes
0 Cap on

	Temperature	O ₂ rate	CO ₂ rate	Preservation time
Classic Cold	+2°C	21 %	0 %	15 d
Janny MT CA	+2°C	1 to 5 %	1 to 5 %	45 d

Ideal curve

Gas content monitoring for Cardoon



Levels are read using the Tiempo Test measuring instrument

Varieties:

Improved White

Puvis

Vert Vaulx Velin

Alger's Red

Sté JANNY MT SAS

La Condemine 71260 Péronne France

Tel : 0033.3.85.23.96.20

Fax : 0033.3.85.36.96.58 / 0033.3.85.23.96.29

Mail : jannymt@jannymt.com

Web : www.jannymt.com