



JANNY MT



Flexible Installations of Controlled Atmosphere Units

Chinese leave



Varieties:

Pé-Tsaï

...

Chinese leave


- Storing during production peaks
- Maintaining flavour quality
- Maintaining colour and firmness
- Avoiding weight loss (less than 1% after 90 d)
- Minimize losses during long storage
- Limiting rot, black spots and preventing excessive leaf peeling when taken out of extended preservation

Picking recommendations

- Match maturity to the desired preservation time
- Pick extremely healthy products
- Do not store products after rain
- Favour small calibres

Key figures

1. Use Janny MT module 610 - 3
2. Fill module with **150 kg loose**
3. Put in cold storage at **+0°C**
4. **Wait until the products are cooled to heart**
5. Close the lid
6. Check gas stabilisation
7. Possible preservation up to **90 days**
8. Open modules
9. Let the module **12 hours with open lid in cold room for rehabilitation**
10. **Use Sorbacal in order to minimise CO₂ rate**




3 Open membranes
3 Caps on

	Temperature	O ₂ rate	CO ₂ rate	Preservation time
Classic Cold	+0°C	21 %	0 %	60 d
Janny MT CA	+0°C	1 to 5 %	1 to 2 %	90 d

Ideal curve

Gas content monitoring for Chinese leave



Levels are read using the Tiempo Test measuring instrument

Sté JANNY MT SAS
 La Condemine 71260 Péronne France
 Tel : 0033.3.85.23.96.20
 Fax : 0033.3.85.36.96.58 / 0033.3.85.23.96.29
 Mail : jannymt@jannymt.com
 Web : www.jannymt.com