# Leek





## Varieties:

Sté JANNY MT SAS La Condemine 71260 Péronne France Tel : 0033.3.85.23.96.20 Fax : 0033.3.85.36.96.58 / 0033.3.85.23.96.29

> Mail: jannymt@jannymt.com Web: www.jannymt.com



### « Janny MT preservation »

### Leek


- Storing during production peaks
- Maintaining the sanitary and gustative quality
- Maintaining colour and firmness
- Avoiding weight loss (less than 1% after 45 d)
- Store before frost
- Preserving before bolting when the plants start growing again

#### Picking recommendations

- Match maturity to the desired preservation time
- Pick extremely healthy products
- Do not store products after rain
- Store washed leeks, with healthy leaves and with a net cut

#### Key figures

- 1. Use Janny MT module 610 6
- 2. Fill module with **150 kg loose**
- 3. Put in cold storage at +0°C
- 4. Wait until the products are cooled to heart
- 5. Close the lid
- 6. Check gas stabilisation
- 7. Possible preservation up to 45 days
- 8. Open modules
- 9. Let the module 12 hours with open lid in cold room for rehabilitation




6 Open membranes

0 Cap on

	Temperature	O <sub>2</sub> rate	CO <sub>2</sub> rate	Preservation time
Classic Cold	+0°C	21 %	0 %	15 to 21 d
Janny MT CA	+0°C	1 to 5 %	1 to 5 %	45 d

#### Ideal curve

# Gas content monitoring for Leek







Levels are read using the Tiempo Test measuring instrument