



JANNY MT



Flexible Installations of Controlled Atmosphere Units

« Janny MT preservation »

Early potato



Early potato


- Lengthen sales period
- Avoiding weight loss (less than 2% after 90 d)
- Maintaining the sanitary and gustative quality
- Maintaining colour and freshness
- Minimize losses during long storage

Picking recommendations

- Match maturity to the desired preservation time
- Pick extremely healthy products
- Do not store products after rain
- Favour small calibres

Key figures

1. Use Janny MT module 610 - 6
2. Fill module with **300 kg loose or 30 trays**
3. Put in cold storage at **+3.5°C**
4. **Wait until the products are cooled to heart**
5. Close the lid
6. Check gas stabilisation
7. Possible preservation up to **90 days**
8. Open modules
9. Let the module **12 hours with open lid in cold room for rehabilitation**




6 Open membranes
0 Cap on

	Temperature	O ₂ rate	CO ₂ rate	Preservation time
Classic Cold	+3.5°C	21 %	0 %	15 to 21 d
Janny MT CA	+3.5°C	1 to 5 %	1 to 5 %	90 d

Ideal curve

Gas content monitoring for Early potato



Levels are read using the Tiempo Test measuring instrument

Varieties:
Délicatesse
Rate du Touquet

...

Sté JANNY MT SAS
La Condemine 71260 Péronne France
Tel : 0033.3.85.23.96.20
Fax : 0033.3.85.36.96.58 / 0033.3.85.23.96.29
Mail : jannymt@jannymt.com
Web : www.jannymt.com